

RAYSTATION 8B

DICOM Conformance Statement Accuray Driver



1 OVERVIEW

This document specifies the DICOM interface for the treatment management system (TMS) Ray Treatment Tomo driver with Accuray TomoTherapy and Radixact Treatment Delivery Devices (TDD).

1.1 NETWORK SERVICES

| SOP Class Name | SOP Class UID | Provider of Service (SCP) | User of Service (SCU) |
|---|-------------------------------|---------------------------|-----------------------|
| Transfer | | | |
| CT Image Storage | 1.2.840.10008.5.1.4.1.1.2 | Yes | No |
| RT Beams Treatment Record Storage | 1.2.840.10008.5.1.4.1.1.481.4 | Yes | No |
| Spatial Registration (REG) Storage | 1.2.840.10008.5.1.4.1.1.66.1 | Yes | No |
| Query/Retrieve | | | |
| Study Root Query/Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Yes | No |
| Workflow Management | | | |
| Unified Procedure Step - Push SOP Class-Trial (Retired) | 1.2.840.10008.5.1.4.34.4.1 | Yes | No |
| Unified Procedure Step - Pull SOP Class-Trial (Retired) | 1.2.840.10008.5.1.4.34.4.3 | Yes | No |
| Verification | | | |
| Verification SOP Class | 1.2.840.10008.1.1 | Yes | No |

3 INTRODUCTION

3.1 REVISION HISTORY

| Date | Version | Comment |
|------------|---------|---|
| 2018-12-13 | 1.0 | Tomo Driver DCS for RayStation Release 8B |
| 2019-02-14 | 2.0 | Tomo Driver DCS for RayStation Release 8B SP1 |

3.2 AUDIENCE

This document is written for users that need to understand how Tomo will integrate into their healthcare facility. This includes both those responsible for overall imaging network policy and architecture, as well as integrators who need to have a detailed understanding of the DICOM features of the product. This document contains some basic DICOM definitions so that any reader may understand how this product implements DICOM features. However, integrators are expected to fully understand all the DICOM terminology, how the tables in this document relate to the product's functionality, and how that functionality integrates with other devices that support compatible DICOM features.

3.3 REMARKS

This document is written for users that need to understand how RayTreatment Driver Tomo will integrate into their healthcare facility. This includes both those responsible for overall imaging network policy and architecture, as well as integrators who need to have a detailed understanding of the DICOM features of the product. This document contains some basic DICOM definitions so that any reader may understand how this product implements DICOM features. However, integrators are expected to fully understand all the DICOM terminology, how the tables in this document relate to the product's functionality, and how that functionality integrates with other devices that support compatible DICOM features.

This Conformance Statement is not supposed to replace validation with other DICOM equipment to ensure proper exchange of intended information. In fact, the user should be aware of the following important issues:

- The comparison of different Conformance Statements is just the first step towards assessing interconnectivity and interoperability between the product and other DICOM conformant equipment.
- Test procedures should be defined and executed to validate the required level of interoperability with specific compatible DICOM equipment, as established by the healthcare facility.

3.3.1 Interoperability validation needed

When using RayTreatment Driver Tomo together with other software, the DICOM conformance statements must be compared and relevant validation tests run. The DICOM standard by itself does not guarantee interoperability. The Conformance Statement does, however, facilitate a first-level comparison for interoperability between different applications supporting compatible DICOM functionality. RaySearch is also active within the IHE-RO. Contact RaySearch for more info regarding adherence to IHE-RO profiles.

3.3.2 DICOM revision

The module tables listed in the last two chapters are based on part 3 of the DICOM-standard revision 2009. For extra clarity all attributes in the referenced modules have been listed, even the ones that are not used by Tomo.

3.4 TERMS AND DEFINITIONS

Informal definitions are provided for the following terms used in this Conformance Statement. The DICOM Standard is the authoritative source for formal definitions of these terms.

Abstract Syntax – the information agreed to be exchanged between applications, generally equivalent to a Service/Object Pair (SOP) Class. Examples: Verification SOP Class, Modality Worklist Information Model Find SOP Class, Computed Radiography Image Storage SOP Class.

Application Entity (AE) – an end point of a DICOM information exchange, including the DICOM network or media interface software; i.e., the software that sends or receives DICOM information objects or messages. A single device may have multiple Application Entities.

Application Entity Title – the externally known name of an Application Entity, used to identify a DICOM application to other DICOM applications on the network.

Application Context – the specification of the type of communication used between Application Entities. Example: DICOM network protocol.

Association – a network communication channel set up between Application Entities.

Attribute – a unit of information in an object definition; a data element identified by a tag. The information may be a complex data structure (Sequence), itself composed of lower level data elements. Examples: Patient ID (0010,0020), Accession Number (0008,0050), Photometric Interpretation (0028,0004), Procedure Code Sequence (0008,1032).

Information Object Definition (IOD) – the specified set of Attributes that comprise a type of data object; does not represent a specific instance of the data object, but rather a class of similar data objects that have the same properties. The Attributes may be specified as Mandatory (Type 1), Required but possibly unknown (Type 2),

or Optional (Type 3), and there may be conditions associated with the use of an Attribute (Types 1C and 2C). Examples: MR Image IOD, CT Image IOD, Print Job IOD.

Joint Photographic Experts Group (JPEG) – a set of standardized image compression techniques, available for use by DICOM applications.

Module – a set of Attributes within an Information Object Definition that are logically related to each other. Example: Patient Module includes Patient Name, Patient ID, Patient Birth Date, and Patient Sex.

Negotiation – first phase of Association establishment that allows Application Entities to agree on the types of data to be exchanged and how that data will be encoded.

Protocol Data Unit (PDU) – a packet (piece) of a DICOM message sent across the network. Devices must specify the maximum size packet they can receive for DICOM messages.

Service Class Provider (SCP) – role of an Application Entity that provides a DICOM network service; typically, a server that performs operations requested by another Application Entity (Service Class User). Examples: Picture Archiving and Communication System (image storage SCP, and image query/retrieve SCP), Radiology Information System (modality worklist SCP).

Service Class User (SCU) – role of an Application Entity that uses a DICOM network service; typically, a client. Examples: imaging modality (image storage SCU, and modality worklist SCU), imaging workstation (image query/retrieve SCU)

Service/Object Pair (SOP) Class – the specification of the network or media transfer (service) of a particular type of data (object); the fundamental unit of DICOM interoperability specification. Examples: Ultrasound Image Storage Service, Basic Grayscale Print Management.

Service/Object Pair (SOP) Instance – an information object; a specific occurrence of information exchanged in a SOP Class. Examples: a specific x-ray image.

Tag – a 32-bit identifier for a data element, represented as a pair of four digit hexadecimal numbers, the “group” and the “element”. If the “group” number is odd, the tag is for a private (manufacturer-specific) data element. Examples: (0010,0020) [Patient ID], (07FE,0010) [Pixel Data], (0019,0210) [private data element]

Transfer Syntax – the encoding used for exchange of DICOM information objects and messages. Examples: JPEG compressed (images), little endian explicit value representation.

Unique Identifier (UID) – a globally unique “dotted decimal” string that identifies a specific object or a class of objects; an ISO-8824 Object Identifier. Examples: Study Instance UID, SOP Class UID, SOP Instance UID.

Value Representation (VR) – the format type of an individual DICOM data element, such as text, an integer, a person’s name, or a code. DICOM information objects can be transmitted with either explicit identification of the type of each data element (Explicit VR), or without explicit identification (Implicit VR); with Implicit VR, the receiving application must use a DICOM data dictionary to look up the format of each data element.

3.5 BASICS OF DICOM COMMUNICATION

This section describes terminology used in this Conformance Statement for the non-specialist. This section is not a substitute for training about DICOM, and it makes many simplifications about the meanings of DICOM terms.

Two Application Entities (devices) that want to communicate with each other over a network using DICOM protocol must first agree on several things during an initial network “handshake”. One of the two devices must initiate an Association (a connection to the other device), and ask if specific services, information, and encoding can be supported by the other device (Negotiation).

DICOM specifies a number of network services and types of information objects, each of which is called an Abstract Syntax for the Negotiation. DICOM also specifies a variety of methods for encoding data, denoted Transfer Syntaxes. The Negotiation allows the initiating Application Entity to propose combinations of Abstract Syntax and Transfer Syntax to be used on the Association; these combinations are called Presentation Contexts. The receiving Application Entity accepts the Presentation Contexts it supports.

For each Presentation Context, the Association Negotiation also allows the devices to agree on Roles – which one is the Service Class User (SCU - client) and which is the Service Class Provider (SCP - server). Normally the device initiating the connection is the SCU, i.e., the client system calls the server, but not always.

The Association Negotiation finally enables exchange of maximum network packet (PDU) size, security information, and network service options (called Extended Negotiation information). The Application Entities, having negotiated the Association parameters, may now commence exchanging data. Common data exchanges include queries for worklists and lists of stored images, transfer of image objects and analyses (structured reports), and sending images to film printers. Each exchangeable unit of data is formatted by the sender in accordance with the appropriate Information Object Definition, and sent using the negotiated Transfer Syntax. There is a Default Transfer Syntax that all systems must accept, but it may not be the most efficient for some use cases. Each transfer is explicitly acknowledged by the receiver with a Response Status indicating success, failure, or that query or retrieve operations are still in process.

Two Application Entities may also communicate with each other by exchanging media (such as a CD-R). Since there is no Association Negotiation possible, they both use a Media Application Profile that specifies “pre-negotiated” exchange media format, Abstract Syntax, and Transfer Syntax.

3.6 ABBREVIATIONS

| Name | Meaning |
|--------------|--|
| AE | Application Entity |
| CT | Computed Tomography |
| DICOM | Digital Imaging and Communications in Medicine |
| IHE / IHE-RO | Integrating the Healthcare Enterprise. IHE-RO deals with integrating Radiation Oncology. |
| IOD | Information Object Definition |
| JPEG | Joint Photographic Experts Group |
| MR | Magnetic Resonance Imaging |
| PACS | Picture Archiving and Communication System |
| PET | Positron Emission Tomography |
| PTS | Proton Planning System (used by IBA) |
| RT | Radiotherapy |
| SCP | Service Class Provider |
| SCU | Service Class User |
| SOP | Service-Object Pair |
| TDD | Treatment Delivery Device |
| TMS | Treatment Management System |
| TPS | Treatment Planning System |

3.7 REFERENCES

- NEMA PS3 Digital Imaging and Communications in Medicine (DICOM) Standard, available free at <http://medical.nema.org/>


4 NETWORKING

The following diagram illustrates the application dataflow between RayTreatment Driver Tomo and the specific TDD.

4.1 IMPLEMENTATION MODEL

4.1.1 Application data flow

The following diagram illustrates the application data flow between RayTreatment Driver Tomo and the specific TDD.



The scenario starts with a C-FIND query for Unified Procedure Steps from the remote client. The client can then take responsibility for the UPS by setting it to IN PROGRESS. Once the UPS is IN PROGRESS the RayTreatment Driver Tomo only allows requests corresponding to the current session until the session has been completed in the application.

4.1.2 Functional Definition of AEs

4.1.2.1 Functional Definition of "Tomo Application Entity"

The following operations are supported:

CT Image

- C-STORE for setup CT images
- C-MOVE for planning CT images

Spartial Registration (SR0)

- C-STORE for registration between setup and planning images.

RT Image

- C-STORE for setup RT Images.

RT Structure Set

- C-MOVE for planning RT Structure Set.

RT Beams Treatment Record

- C-STORE for delivery result.
- C-MOVE for previous delivered results.

Beams Delivery Instructions

- C-MOVE for BDIs related to the Unified Procedure Step.

Modality Performed Procedure Step - PULL

- C-FIND for worklist query.
- N-ACTION for UPS status changes.
- N-SET for progress update.

Verification

- C-ECHO for connection verification

4.1.3 Sequence of Real World Activities

4.1.3.1 Prepare session

Once the patient is checked in to the session, Unified Procedure Steps will be created and available for Worklist queries.

4.1.3.2 Manual cancellation

The procedure step can be canceled by the user in the application. Further requests relation to the session will be rejected.

4.1.3.3 Complete session

All sessions, including canceled sessions, needs to be completed by the user in the application before another session can be started.

4.2 AE SPECIFICATIONS:

4.2.1 Raytreatment Driver Tomo Application Entity

4.2.1.1 SOP Classes

| SOP Class Name | SOP Class UID | Provider of Service (SCP) | User of Service (SCU) |
|---|-------------------------------|---------------------------|-----------------------|
| Transfer | | | |
| CT Image Storage | 1.2.840.10008.5.1.4.1.1.2 | Yes | No |
| RT BeamsTreatment Record Storage | 1.2.840.10008.5.1.4.1.1.481.4 | Yes | No |
| Spatial Registration (REG) Storage | 1.2.840.10008.5.1.4.1.1.66.1 | Yes | No |
| Query/Retrieve | | | |
| Study Root Query/Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Yes | No |
| Workflow Management | | | |
| Unified Procedure Step - Push SOP Class-Trial (Retired) | 1.2.840.10008.5.1.4.34.4.1 | Yes | No |
| Unified Procedure Step - Pull SOP Class-Trial (Retired) | 1.2.840.10008.5.1.4.34.4.3 | Yes | No |
| Verification | | | |
| Verification SOP Class | 1.2.840.10008.1.1 | Yes | No |

4.2.1.2 Association Policies

Not applicable

4.2.1.3 General

The DICOM standard Application context shall be specified.

| | |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

The maximum PDU size in not configurable and is set to 16384 for SCU and unlimited for SCP.

4.2.1.4 Number of Associations

Any number of incoming concurrent associations are accepted.

4.2.1.5 Asynchronous Nature

RayTreatment Driver Tomo does not support asynchronous communication (multiple outstanding transactions over a single Association).

4.2.1.6 Implementation Identity Information

Not applicable

4.2.1.7 Association Initiation Policy

The implementation for this Application Entity is:

| | |
|-----------------------------|---------------------|
| Implementation Class UID | 1.3.6.1.4.1.30071.8 |
| Implementation Version Name | fo-dicom 3.0.5 |

4.2.1.8 Activity C-ECHO

4.2.1.8.1 Description and Sequencing of Activities

A C-ECHO request can always be sent to the Tomo driver.

4.2.1.8.2 Accepted Presentation Context

| Abstract Syntax Name | Abstract Syntax UID | Transfer Syntax name | Transfer Syntax UID | Role | Extended Negotiation |
|----------------------|---------------------|---------------------------|---------------------|------|----------------------|
| Verification | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

4.2.1.9 Activity C-FIND

4.2.1.9.1 Description and Sequencing of Activities

Not applicable

4.2.1.9.2 Accepted Presentation Context

| Abstract Syntax Name | Abstract Syntax UID | Transfer Syntax name | Transfer Syntax UID | Role | Extended Negotiation |
|---|----------------------------|---------------------------|---------------------|------|----------------------|
| Unified Procedure Step - Pull SOP Class - Trial (Retired) | 1.2.840.10008.5.1.4.34.4.3 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

4.2.1.10 Activity C-MOVE

4.2.1.10.1 Description and Sequencing of Activities

Not applicable

4.2.1.10.2 Accepted Presentation Context

RayTreatment Driver Tomo does not support C-MOVE requests.

4.2.1.11 Activity C-STORE

4.2.1.11.1 Description and Sequencing of Activities

Not applicable

4.2.1.11.2 Accepted Presentation Context

| Abstract Syntax Name | Abstract Syntax UID | Transfer Syntax name | Transfer Syntax UID | Role | Extended Negotiation |
|--------------------------|-------------------------------|---------------------------|---------------------|------|----------------------|
| CT Image Storage | 1.2.840.10008.5.1.4.1.1.2 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| RT Structure Set Storage | 1.2.840.10008.5.1.4.1.1.481.3 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |

| | | | | | |
|-----------------------------------|-------------------------------|---------------------------|---------------------|-----|------|
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| Spatial Registration Storage | 1.2.840.10008.5.1.4.1.1.66.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| RT Beams Treatment Record Storage | 1.2.840.10008.5.1.4.1.1.481.4 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

4.2.1.11.3 Status Response

| Service Status | Further meaning | Error Code | Reason |
|----------------|---------------------------|------------|---|
| Failure | Storage Cannot Understand | Cxxx | Cannot find session or validation failed. |
| | SOP class not supported | 0122 | SOP class not supported. |
| Success | Success | 0000 | |

4.2.1.12 Activity N-ACTION

4.2.1.12.1 Description and Sequencing of Activities

Not applicable

4.2.1.12.2 Accepted Presentation Context

| Abstract Syntax Name | Abstract Syntax UID | Transfer Syntax name | Transfer Syntax UID | Role | Extended Negotiation |
|---|----------------------------|---------------------------|---------------------|------|----------------------|
| Unified Procedure Step - Pull SOP Class - Trial (Retired) | 1.2.840.10008.5.1.4.34.4.3 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

4.2.1.12.3 Status Response

| Service Status | Further meaning | Error Code | Reason |
|----------------|------------------------------|------------|---|
| Refused | QueryRetrieveUnableToProcess | C000 | Unknown SOP Instance UID. |
| | NoLongerUpdateUps | C300 | The UPS may no longer be updated. |
| | IncorrectTransactionUid | C301 | The correct Transaction UID was not provided. |
| | AlreadyInProgress | C302 | The UPS is already IN PROGRESS |
| | IncorrectTransactionUid | C301 | The correct Transaction UID was not provided. |
| | SopInstanceUidDoesNotExists | C307 | Specified SOP Instance UID does not exist or is nota UPS Instance managed by this SCP |
| Failure | AlreadyCanceled | B304 | The UPS is already in the requested state of CANCELED. |
| | AlreadyCompleted | B306 | The UPS is already in the requested state of COMPLETED. |
| Success | Success | 0000 | |

4.2.1.13 Activity N-SET

4.2.1.13.1 Description and Sequencing of Activities

Not applicable

4.2.1.13.2 Accepted Presentation Context

| Abstract Syntax Name | Abstract Syntax UID | Transfer Syntax name | Transfer Syntax UID | Role | Extended Negotiation |
|---|----------------------------|---------------------------|---------------------|------|----------------------|
| Unified Procedure Step - Pull SOP Class - Trial (Retired) | 1.2.840.10008.5.1.4.34.4.3 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |

| | | | | | |
|--|--|---------------------------|---------------------|--|--|
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

4.2.1.13.3 Status Response

| Service Status | Further meaning | Error Code | Reason |
|----------------|------------------------------|------------|---|
| Refused | QueryRetrieveUnableToProcess | C000 | Unknown SOP Instance UID. |
| | NoLongerUpdateUps | C300 | The UPS may no longer be updated. |
| | IncorrectTransactionUid | C301 | The correct Transaction UID was not provided. |
| | SopInstanceUidDoesNotExists | C307 | Specified SOP Instance UID does not exist or is nota UPS Instance managed by this SCP |
| Success | Success | 0000 | |

5 MEDIA INTERCHANGE

Not applicable

6 TRANSFORMATION OF DICOM TO CDA

Not applicable

7 SUPPORT OF CHARACTER SETS

RayTreatment Driver Tomo support the following charactersets in addition to the default

- ISO_IR 192

8 SECURITY

8.1 SECURITY PROFILES

No Security Profiles are supported.

8.2 ASSOCIATION LEVEL SECURITY

RayTreatment Driver Tomo checks the following values for validation of received Association Open Requests:

- Called AE Title.

8.3 APPLICATION LEVEL SECURITY

None supported.

9 ANNEXES

9.1 IOD CONTENTS

9.1.1 Created SOP Instance(s)

9.1.1.1 RT Plan IOD

| IE | Module | Used |
|-----------|-------------------------------------|------|
| Patient | Patient Module | No |
| Study | General Study Module | No |
| Series | RT Series Module | No |
| Equipment | General Equipment Module | No |
| Plan | RT General Plan Module | No |
| | RT Beams Module | No |
| | RT Brachy Application Setups Module | No |
| | SOP Common Module | No |

9.1.1.2 RT Beams Treatment Record IOD

| IE | Module | Used |
|------------------|------------------------------------|------|
| Patient | Patient Module | Yes |
| Study | General Study Module | Yes |
| Series | RT Series Module | Yes |
| Equipment | General Equipment Module | Yes |
| Treatment Record | RT General Treatment Record Module | Yes |
| | RT Patient Setup Module | Yes |
| | RT Treatment Machine Record Module | Yes |
| | RT Beams Session Record Module | Yes |
| | SOP Common Module | Yes |

9.1.1.2.1 Patient Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------------|-------------|----|------|---------|
| Patient's Name | (0010,0010) | PN | 2 | |
| Patient ID | (0010,0020) | LO | 2 | |
| Patient's Birth Date | (0010,0030) | DA | 2 | |
| Patient's Sex | (0010,0040) | CS | 2 | |

9.1.1.2.2 General Study Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------------------|-------------|----|------|---------|
| Study Instance UID | (0020,0000) | UI | 1 | |
| Study Date | (0008,0020) | DA | 2 | |
| Study Time | (0008,0030) | TM | 2 | |
| Referring Physician's Name | (0008,0090) | PN | 2 | |
| Study ID | (0020,0010) | SH | 2 | |
| Study Description | (0008,1030) | LO | 3 | |

9.1.1.2.3 RT Series Module

| Attribute name | Tag | Vr | Type | Comment |
|---------------------|-------------|----|------|------------------|
| Modality | (0008,0060) | CS | 1 | Always RTRECORD. |
| Series Instance UID | (0020,000E) | UI | 1 | |
| Series Number | (0020,0011) | IS | 2 | |
| Operators' Name | (0008,1070) | PN | 2 | |

9.1.1.2.4 General Equipment Module

| Attribute name | Tag | Vr | Type | Comment |
|---------------------------|-------------|----|------|---------|
| Manufacturer | (0008,0070) | LO | 2 | |
| Manufacturer's Model Name | (0008,1090) | LO | 3 | |
| Software Version(s) | (0018,1020) | LO | 3 | |

9.1.1.2.5 RT General Treatment Record Module

| Attribute name | Tag | Vr | Type | Comment |
|--------------------------------------|-------------|----|------|---------|
| Instance Number | (0020,0013) | IS | 1 | |
| Treatment Date | (3008,0250) | DA | 2 | |
| Treatment Time | (3008,0251) | TM | 2 | |
| Referenced RT Plan Sequence | (300C,0002) | SQ | 2 | |
| >Referenced SOP Class UID | (0008,1150) | UI | 1 | |
| >Referenced SOP Instance UID | (0008,1155) | UI | 1 | |
| Referenced Treatment Record Sequence | (3008,0030) | SQ | 3 | |
| >Referenced SOP Class UID | (0008,1150) | UI | 1 | |
| >Referenced SOP Instance UID | (0008,1155) | UI | 1 | |

9.1.1.2.6 RT Patient Setup Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------------|-------------|----|------|---------|
| Patient Setup Sequence | (300A,0180) | SQ | 1 | |
| >Patient Setup Number | (300A,0182) | IS | 1 | |
| >Patient Position | (0018,5100) | CS | 1C | |
| >Setup Technique | (300A,01B0) | CS | 3 | |

9.1.1.2.7 RT Treatment Machine Record Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------------------|-------------|----|------|---------|
| Treatment Machine Sequence | (300A,0206) | SQ | 1 | |
| >Manufacturer | (0008,0070) | LO | 2 | |
| >Institution Name | (0008,0080) | LO | 2 | |
| >Manufacturer's Model Name | (0008,1090) | LO | 2 | |
| >Device Serial Number | (0018,1000) | LO | 2 | |

9.1.1.2.8 RT Beams Session Record Module

| Attribute name | Tag | Vr | Type | Comment |
|---------------------------------|-------------|----|------|---------|
| Number of Fractions Planned | (300A,0078) | IS | 2 | |
| Primary Dosimeter Unit | (300A,00B3) | CS | 1 | |
| Treatment Session Beam Sequence | (3008,0020) | SQ | 1 | |
| >Referenced Beam Number | (300C,0006) | IS | 3 | |
| >Beam Name | (300A,00C2) | LO | 3 | |

| | | | | |
|---|-------------|----|----|--|
| >Beam Type | (300A,00C4) | CS | 1 | |
| >Radiation Type | (300A,00C6) | CS | 1 | |
| >Beam Limiting Device Leaf Pairs Sequence | (3008,00A0) | SQ | 1 | |
| >RT Beam Limiting Device Type | (300A,00B8) | CS | 1 | |
| >Number of Leaf/Jaw Pairs | (300A,00BC) | IS | 1 | |
| >Referenced Patient Setup Number | (300C,006A) | IS | 3 | |
| >Number of Wedges | (300A,00D0) | IS | 1 | |
| >Number of Compensators | (300A,00E0) | IS | 2 | |
| >Number of Boli | (300A,00ED) | IS | 2 | |
| >Number of Blocks | (300A,00F0) | IS | 2 | |
| >Current Fraction Number | (3008,0022) | IS | 2 | |
| >Treatment Delivery Type | (300A,00CE) | CS | 2 | |
| >Treatment Termination Status | (3008,002A) | CS | 1 | |
| >Treatment Verification Status | (3008,002C) | CS | 2 | |
| >Specified Primary Meterset | (3008,0032) | DS | 3 | |
| >Delivered Primary Meterset | (3008,0036) | DS | 3 | |
| >Number of Control Points | (300A,0110) | IS | 1 | |
| >Control Point Delivery Sequence | (3008,0040) | SQ | 1 | |
| >Referenced Control Point Index | (300C,00F0) | IS | 3 | |
| >Treatment Control Point Date | (3008,0024) | DA | 1 | |
| >Treatment Control Point Time | (3008,0025) | TM | 1 | |
| >Specified Meterset | (3008,0042) | DS | 2 | |
| >Delivered Meterset | (3008,0044) | DS | 1 | |
| >Nominal Beam Energy | (300A,0114) | DS | 3 | |
| >Wedge Position Sequence | (300A,0116) | SQ | 3 | |
| >Referenced Wedge Number | (300C,00C0) | IS | 1 | |
| >Wedge Position | (300A,0118) | CS | 1 | |
| >Gantry Angle | (300A,011E) | DS | 1C | |
| >Gantry Rotation Direction | (300A,011F) | CS | 1C | |
| >Beam Limiting Device Angle | (300A,0120) | DS | 1C | |
| >Beam Limiting Device Rotation Direction | (300A,0121) | CS | 1C | |
| >Patient Support Angle | (300A,0122) | DS | 1C | |
| >Patient Support Rotation Direction | (300A,0123) | CS | 1C | |
| >Table Top Pitch Angle | (300A,0140) | FL | 1C | |
| >Table Top Pitch Rotation Direction | (300A,0142) | CS | 1C | |
| >Table Top Roll Angle | (300A,0144) | FL | 1C | |
| >Table Top Roll Rotation Direction | (300A,0146) | CS | 1C | |
| >Table Top Vertical Position | (300A,0128) | DS | 2C | |
| >Table Top Longitudinal Position | (300A,0129) | DS | 2C | |
| >Table Top Lateral Position | (300A,012A) | DS | 2C | |

9.1.1.2.9 SOP Common Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------|-------------|----|------|---------|
| SOP Class UID | (0008,0016) | UI | 1 | |

| | | | | |
|------------------|-------------|----|---|--|
| SOP Instance UID | (0008,0018) | UI | 1 | |
|------------------|-------------|----|---|--|

9.1.1.3 RT Beams Delivery Instruction RETIRED IOD

| IE | Module | Used |
|-----------|--------------------------------------|------|
| Patient | Patient Module | Yes |
| | Clinical Trial Subject Module | No |
| Study | General Study Module | Yes |
| | Patient Study Module | No |
| | Clinical Trial Study Module | No |
| Series | General Series Module | Yes |
| | Clinical Trial Series Module | No |
| Equipment | General Equipment Module | Yes |
| Plan | RT Beams Delivery Instruction Module | Yes |
| | Common Instance Reference Module | No |
| | General Reference Module | No |
| | SOP Common Module | Yes |

9.1.1.3.1 Patient Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------------|-------------|----|------|---------------------------|
| Patient's Name | (0010,0010) | PN | 2 | |
| Patient ID | (0010,0020) | LO | 2 | |
| Patient's Birth Date | (0010,0030) | DA | 2 | |
| Patient's Sex | (0010,0040) | CS | 2 | Possible values: M, F, O. |

9.1.1.3.2 General Study Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------------------|-------------|----|------|---------|
| Study Instance UID | (0020,0000) | UI | 1 | |
| Study Date | (0008,0020) | DA | 2 | |
| Study Time | (0008,0030) | TM | 2 | |
| Referring Physician's Name | (0008,0090) | PN | 2 | |
| Study ID | (0020,0010) | SH | 2 | |
| Accession Number | (0008,0050) | SH | 2 | |
| Study Description | (0008,1030) | LO | 3 | |

9.1.1.3.3 General Series Module

| Attribute name | Tag | Vr | Type | Comment |
|---------------------|-------------|----|------|--------------|
| Modality | (0008,0060) | CS | 1 | Always PLAN. |
| Series Instance UID | (0020,000E) | UI | 1 | |
| Series Number | (0020,0011) | IS | 2 | |
| Series Date | (0008,0021) | DA | 3 | |
| Series Time | (0008,0031) | TM | 3 | |
| Series Description | (0008,103E) | LO | 3 | |

9.1.1.3.4 General Equipment Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------|-------------|----|------|---------|
| Manufacturer | (0008,0070) | LO | 2 | |
| Station Name | (0008,1010) | SH | 3 | |

9.1.1.3.5 RT Beams Delivery Instruction Module

| Attribute name | Tag | Vr | Type | Comment |
|--|-------------|----|------|--|
| Beam Task Sequence | (0074,1020) | SQ | 1 | |
| >Beam Task Type | (0074,1022) | CS | 1 | Always <ul style="list-style-type: none"> TREAT - Treat |
| >Treatment Delivery Type | (300A,00CE) | CS | 1 | Possible values: TREATMENT, CONTINUATION. |
| >Continuation Start Meterset | (0074,0120) | FD | 1C | |
| >Continuation End Meterset | (0074,0121) | FD | 1C | |
| >Current Fraction Number | (3008,0022) | IS | 1 | |
| >Referenced Beam Number | (300C,0006) | IS | 1 | |
| >Table Top Vertical Setup Displacement | (300A,01D2) | DS | 2 | |
| >Table Top Longitudinal Setup Displacement | (300A,01D4) | DS | 2 | |
| >Table Top Lateral Setup Displacement | (300A,01D6) | DS | 2 | |
| >Referenced RT Plan Sequence | (300C,0002) | SQ | 3 | |
| >Referenced Series Sequence | (0008,1115) | SQ | 1C | |
| >Series Instance UID | (0020,000E) | UI | 1 | |
| >Retrieve AE Title | (0008,0054) | AE | | |
| >Referenced SOP Sequence | (0008,1199) | SQ | | |
| >Referenced SOP Class UID | (0008,1150) | UI | 1 | |
| >Referenced SOP Instance UID | (0008,1155) | UI | 1 | |
| >Study Instance UID | (0020,000D) | UI | 1 | |
| Omitted Beam Task Sequence | (300C,0111) | SQ | 3 | |
| >Referenced Beam Number | (300C,0006) | IS | 1 | |
| >Reason for Omission | (300C,0112) | CS | 1 | Always ALREADY_TREATED. |

9.1.1.3.6 SOP Common Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------------|-------------|----|------|---------|
| SOP Class UID | (0008,0016) | UI | 1 | |
| SOP Instance UID | (0008,0018) | UI | 1 | |
| Specific Character Set | (0008,0005) | CS | 1C | |
| Instance Creation Date | (0008,0012) | DA | 3 | |
| Instance Creation Time | (0008,0013) | TM | 3 | |

9.1.1.4 Unified Procedure Step RETIRED IOD

| IE | Module | Used |
|------------------------|---|------|
| Unified Procedure Step | SOP Common Module | Yes |
| | Unified Procedure Step Relationship Module | Yes |
| | Unified Procedure Step Scheduled Procedure Information Module | Yes |
| | Unified Procedure Step Progress Information Module | Yes |
| | Unified Procedure Step Performed Procedure Information Module | No |
| | Patient Medical Module | No |
| | Transaction Module | Yes |

9.1.1.4.1 SOP Common Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------|-----|----|------|---------|
|----------------|-----|----|------|---------|

| | | | | |
|--------------------------|-------------|----|---|--|
| SOP Class UID | (0008,0016) | UI | 1 | |
| SOP Instance UID | (0008,0018) | UI | 1 | |
| Timezone Offset From UTC | (0008,0201) | SH | 3 | |

9.1.1.4.2 Unified Procedure Step Relationship Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------------------------|-------------|----|------|--|
| Patient's Name | (0010,0010) | PN | | |
| Patient ID | (0010,0020) | LO | | |
| Patient's Birth Date | (0010,0030) | DA | | |
| Patient's Sex | (0010,0040) | CS | | Possible values: <ul style="list-style-type: none"> • M - Male • F - Female • O - Other |
| Referenced Request Sequence | (0040,A370) | SQ | | |
| >Requested Procedure Code Sequence | (0032,1064) | SQ | | |
| >Code Value | (0008,0100) | SH | 1C | |
| >Coding Scheme Designator | (0008,0102) | SH | 1C | |
| >Coding Scheme Version | (0008,0103) | SH | 1C | |
| >Code Meaning | (0008,0104) | LO | 1 | |
| >Requested Procedure ID | (0040,1001) | SH | | |

9.1.1.4.3 Unified Procedure Step Scheduled Procedure Information Module

| Attribute name | Tag | Vr | Type | Comment |
|--|-------------|----|------|--|
| Scheduled Procedure Step Priority | (0074,1200) | CS | | Possible values: <ul style="list-style-type: none"> • HIGH - High • MEDIUM - Medium • LOW - Low |
| Procedure Step Label | (0074,1204) | LO | | |
| Scheduled Station Name Code Sequence | (0040,4025) | SQ | | |
| >Code Value | (0008,0100) | SH | 1C | |
| >Coding Scheme Designator | (0008,0102) | SH | 1C | |
| >Code Meaning | (0008,0104) | LO | 1 | |
| Scheduled Procedure Step Start DateTime | (0040,4005) | DT | | |
| Expected Completion DateTime | (0040,4011) | DT | | |
| Scheduled Workitem Code Sequence | (0040,4018) | SQ | | |
| >Code Value | (0008,0100) | SH | 1C | |
| >Coding Scheme Designator | (0008,0102) | SH | 1C | |
| >Code Meaning | (0008,0104) | LO | 1 | |
| Scheduled Processing Parameters Sequence | (0074,1210) | SQ | | |
| >Value Type | (0040,A040) | CS | 1 | Always <ul style="list-style-type: none"> • TEXT - Text |
| >Concept Name Code Sequence | (0040,A043) | SQ | 1 | |
| >Code Value | (0008,0100) | SH | 1C | |
| >Coding Scheme Designator | (0008,0102) | SH | 1C | |
| >Code Meaning | (0008,0104) | LO | 1 | |

| | | | | |
|------------------------------|-------------|----|----|---|
| >Text Value | (0040,A160) | UT | 1C | |
| Input Information Sequence | (0040,4021) | SQ | | |
| >Study Instance UID | (0020,000D) | UI | 1 | |
| >Referenced Series Sequence | (0008,1115) | SQ | 1C | |
| >Series Instance UID | (0020,000E) | UI | 1 | |
| >Retrieve AE Title | (0008,0054) | AE | | |
| >Referenced SOP Sequence | (0008,1199) | SQ | | |
| >Referenced SOP Class UID | (0008,1150) | UI | 1 | |
| >Referenced SOP Instance UID | (0008,1155) | UI | 1 | |
| Study Instance UID | (0020,000D) | UI | | |
| Input Availability Flag | (0040,4020) | CS | 1 | Possible values: <ul style="list-style-type: none"> • COMPLETE - Complete • INCOMPLETE - Incomplete |

9.1.1.4.4 Unified Procedure Step Progress Information Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------------|-------------|----|------|---|
| Procedure Step State | (0074,1000) | CS | | Possible values: <ul style="list-style-type: none"> • SCHEDULED - Scheduled • IN PROGRESS - In Progress • CANCELED - Canceled • COMPLETED - Completed |

9.1.1.4.5 Transaction Module

| Attribute name | Tag | Vr | Type | Comment |
|-----------------|-------------|----|------|---------|
| Transaction UID | (0008,1195) | UI | 3 | |

9.1.2 Usage of Attributes From Received IODs

9.1.2.1 CT Image IOD

| IE | Module | Used |
|--------------------|---------------------------|------|
| Patient | Patient Module | Yes |
| Study | General Study Module | Yes |
| Series | General Series Module | Yes |
| Frame of Reference | Frame of Reference Module | Yes |
| Equipment | General Equipment Module | No |
| Image | General Image Module | No |
| | Image Plane Module | No |
| | Image Pixel Module | No |
| | Contrast/Bolus Module | No |
| | CT Image Module | Yes |
| | SOP Common Module | Yes |

9.1.2.1.1 Patient Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------------|-------------|----|------|---------|
| Patient's Name | (0010,0010) | PN | 2 | |
| Patient ID | (0010,0020) | LO | 2 | |
| Patient's Birth Date | (0010,0030) | DA | 2 | |

| | | | | |
|---------------|-------------|----|---|--|
| Patient's Sex | (0010,0040) | CS | 2 | |
|---------------|-------------|----|---|--|

9.1.2.1.2 General Study Module

| Attribute name | Tag | Vr | Type | Comment |
|--------------------|-------------|----|------|---------|
| Study Instance UID | (0020,000D) | UI | 1 | |

9.1.2.1.3 General Series Module

| Attribute name | Tag | Vr | Type | Comment |
|---------------------|-------------|----|------|---------|
| Modality | (0008,0060) | CS | 1 | |
| Series Instance UID | (0020,000E) | UI | 1 | |
| Series Date | (0008,0021) | DA | 3 | |
| Series Time | (0008,0031) | TM | 3 | |
| Series Description | (0008,103E) | LO | 3 | |
| Patient Position | (0018,5100) | CS | 2C | |

9.1.2.1.4 Frame of Reference Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------------|-------------|----|------|---------|
| Frame of Reference UID | (0020,0052) | UI | 1 | |

9.1.2.1.5 CT Image Module

| Attribute name | Tag | Vr | Type | Comment |
|---------------------------------|-------------|----|------|---------|
| Image Type | (0008,0008) | CS | 1 | |
| Samples per Pixel | (0028,0002) | US | 1 | |
| Photometric Interpretation | (0028,0004) | CS | 1 | |
| Bits Allocated | (0028,0100) | US | 1 | |
| Bits Stored | (0028,0101) | US | 1 | |
| High Bit | (0028,0102) | US | 1 | |
| Rescale Intercept | (0028,1052) | DS | 1 | |
| Rescale Slope | (0028,1053) | DS | 1 | |
| Table Height | (0018,1130) | DS | 3 | |
| Patient Support Angle | (300A,0122) | DS | 3 | |
| Table Top Pitch Angle | (300A,0140) | FL | 3 | |
| Table Top Roll Angle | (300A,0144) | FL | 3 | |
| Table Top Longitudinal Position | (300A,0129) | DS | 3 | |
| Table Top Lateral Position | (300A,012A) | DS | 3 | |

9.1.2.1.6 SOP Common Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------|-------------|----|------|---------|
| SOP Class UID | (0008,0016) | UI | 1 | |
| SOP Instance UID | (0008,0018) | UI | 1 | |

9.1.2.2 RT Structure Set IOD

| IE | Module | Used |
|-----------|--------------------------|------|
| Patient | Patient Module | Yes |
| Study | General Study Module | Yes |
| Series | RT Series Module | Yes |
| Equipment | General Equipment Module | No |

| | | |
|---------------|----------------------------|-----|
| Structure Set | Structure Set Module | Yes |
| | ROI Contour Module | Yes |
| | RT ROI Observations Module | Yes |
| | SOP Common Module | Yes |

9.1.2.2.1 Patient Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------------|-------------|----|------|---------|
| Patient's Name | (0010,0010) | PN | 2 | |
| Patient ID | (0010,0020) | LO | 2 | |
| Patient's Birth Date | (0010,0030) | DA | 2 | |
| Patient's Sex | (0010,0040) | CS | 2 | |

9.1.2.2.2 General Study Module

| Attribute name | Tag | Vr | Type | Comment |
|--------------------|-------------|----|------|---------|
| Study Instance UID | (0020,0000) | UI | 1 | |

9.1.2.2.3 RT Series Module

| Attribute name | Tag | Vr | Type | Comment |
|---------------------|-------------|----|------|---------|
| Modality | (0008,0060) | CS | 1 | |
| Series Instance UID | (0020,000E) | UI | 1 | |

9.1.2.2.4 Structure Set Module

| Attribute name | Tag | Vr | Type | Comment |
|--|-------------|----|------|---------|
| Structure Set Label | (3006,0002) | SH | 1 | |
| Referenced Frame of Reference Sequence | (3006,0010) | SQ | 3 | |
| >Frame of Reference UID | (0020,0052) | UI | 1 | |
| >RT Referenced Study Sequence | (3006,0012) | SQ | 3 | |
| >Referenced SOP Class UID | (0008,1150) | UI | 1 | |
| >Referenced SOP Instance UID | (0008,1155) | UI | 1 | |
| >RT Referenced Series Sequence | (3006,0014) | SQ | 1 | |
| >Series Instance UID | (0020,000E) | UI | 1 | |
| >Contour Image Sequence | (3006,0016) | SQ | 1 | |
| >Referenced SOP Class UID | (0008,1150) | UI | 1 | |
| >Referenced SOP Instance UID | (0008,1155) | UI | 1 | |
| Structure Set ROI Sequence | (3006,0020) | SQ | 1 | |

9.1.2.2.5 ROI Contour Module

| Attribute name | Tag | Vr | Type | Comment |
|---------------------------|-------------|----|------|---------|
| ROI Contour Sequence | (3006,0039) | SQ | 1 | |
| >Referenced ROI Number | (3006,0084) | IS | 1 | |
| >Contour Sequence | (3006,0040) | SQ | 3 | |
| >Contour Geometric Type | (3006,0042) | CS | 1 | |
| >Number of Contour Points | (3006,0046) | IS | 1 | |
| >Contour Data | (3006,0050) | DS | 1 | |

9.1.2.2.6 RT ROI Observations Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------------------|-------------|----|------|---------|
| RT ROI Observations Sequence | (3006,0080) | SQ | 1 | |

| | | | | |
|-----------------------------------|-------------|----|---|--|
| >Observation Number | (3006,0082) | IS | 1 | |
| >Referenced ROI Number | (3006,0084) | IS | 1 | |
| >RT ROI Interpreted Type | (3006,00A4) | CS | 2 | |
| >ROI Physical Properties Sequence | (3006,00B0) | SO | 3 | |
| >ROI Physical Property | (3006,00B2) | CS | 1 | |
| >ROI Physical Property Value | (3006,00B4) | DS | 1 | |

9.1.2.2.7 SOP Common Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------|-------------|----|------|---------|
| SOP Class UID | (0008,0016) | UI | 1 | |
| SOP Instance UID | (0008,0018) | UI | 1 | |

9.1.2.3 RT Plan IOD

| IE | Module | Used |
|--------------------|-------------------------------------|------|
| Patient | Patient Module | Yes |
| Study | General Study Module | Yes |
| Series | RT Series Module | Yes |
| Frame of Reference | Frame of Reference Module | Yes |
| Equipment | General Equipment Module | No |
| Plan | RT General Plan Module | Yes |
| | RT Patient Setup Module | Yes |
| | RT Fraction Scheme Module | Yes |
| | RT Beams Module | Yes |
| | RT Brachy Application Setups Module | No |
| | SOP Common Module | Yes |

9.1.2.3.1 Patient Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------------|-------------|----|------|---------|
| Patient's Name | (0010,0010) | PN | 2 | |
| Patient ID | (0010,0020) | LO | 2 | |
| Patient's Birth Date | (0010,0030) | DA | 2 | |
| Patient's Sex | (0010,0040) | CS | 2 | |

9.1.2.3.2 General Study Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------------------|-------------|----|------|---------|
| Study Instance UID | (0020,0000) | UI | 1 | |
| Study Date | (0008,0020) | DA | 2 | |
| Study Time | (0008,0030) | TM | 2 | |
| Referring Physician's Name | (0008,0090) | PN | 2 | |
| Study ID | (0020,0010) | SH | 2 | |
| Accession Number | (0008,0050) | SH | 2 | |
| Study Description | (0008,1030) | LO | 3 | |

9.1.2.3.3 RT Series Module

| Attribute name | Tag | Vr | Type | Comment |
|---------------------|-------------|----|------|---------|
| Modality | (0008,0060) | CS | 1 | |
| Series Instance UID | (0020,000E) | UI | 1 | |

9.1.2.3.4 Frame of Reference Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------------|-------------|----|------|---------|
| Frame of Reference UID | (0020,0052) | UI | 1 | |

9.1.2.3.5 RT General Plan Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------|-------------|----|------|---------|
| RT Plan Label | (300A,0002) | SH | 1 | |
| RT Plan Geometry | (300A,000C) | CS | 1 | |

9.1.2.3.6 RT Patient Setup Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------------|-------------|----|------|---------|
| Patient Setup Sequence | (300A,0180) | SQ | 1 | |
| >Patient Setup Number | (300A,0182) | IS | 1 | |
| >Patient Position | (0018,5100) | CS | 1C | |
| >Setup Technique | (300A,01B0) | CS | 3 | |

9.1.2.3.7 RT Fraction Scheme Module

| Attribute name | Tag | Vr | Type | Comment |
|--------------------------------------|-------------|----|------|---------|
| Fraction Group Sequence | (300A,0070) | SQ | 1 | |
| >Fraction Group Number | (300A,0071) | IS | 1 | |
| >Number of Fractions Planned | (300A,0078) | IS | 2 | |
| >Number of Beams | (300A,0080) | IS | 1 | |
| >Number of Brachy Application Setups | (300A,00A0) | IS | 1 | |

9.1.2.3.8 RT Beams Module

| Attribute name | Tag | Vr | Type | Comment |
|--|-------------|----|------|---------|
| Beam Sequence | (300A,00B0) | SQ | 1 | |
| >Beam Number | (300A,00C0) | IS | 1 | |
| >Beam Name | (300A,00C2) | LO | 3 | |
| >Beam Type | (300A,00C4) | CS | 1 | |
| >Radiation Type | (300A,00C6) | CS | 2 | |
| >Beam Limiting Device Sequence | (300A,00B6) | SQ | 1 | |
| >Treatment Delivery Type | (300A,00CE) | CS | 3 | |
| >Number of Wedges | (300A,00D0) | IS | 1 | |
| >Number of Compensators | (300A,00E0) | IS | 1 | |
| >Number of Boli | (300A,00ED) | IS | 1 | |
| >Number of Blocks | (300A,00F0) | IS | 1 | |
| >Number of Control Points | (300A,0110) | IS | 1 | |
| >Control Point Sequence | (300A,0111) | SQ | 1 | |
| >Control Point Index | (300A,0112) | IS | 1 | |
| >Nominal Beam Energy | (300A,0114) | DS | 3 | |
| >Gantry Angle | (300A,011E) | DS | 1C | |
| >Gantry Rotation Direction | (300A,011F) | CS | 1C | |
| >Beam Limiting Device Angle | (300A,0120) | DS | 1C | |
| >Beam Limiting Device Rotation Direction | (300A,0121) | CS | 1C | |
| >Patient Support Angle | (300A,0122) | DS | 1C | |
| >Patient Support Rotation Direction | (300A,0123) | CS | 1C | |

| | | | | |
|-------------------------------------|-------------|----|----|--|
| >Table Top Pitch Angle | (300A,0140) | FL | 1C | |
| >Table Top Pitch Rotation Direction | (300A,0142) | CS | 1C | |
| >Table Top Roll Angle | (300A,0144) | FL | 1C | |
| >Table Top Roll Rotation Direction | (300A,0146) | CS | 1C | |
| >Table Top Vertical Position | (300A,0128) | DS | 2C | |
| >Table Top Longitudinal Position | (300A,0129) | DS | 2C | |
| >Table Top Lateral Position | (300A,012A) | DS | 2C | |
| >Isocenter Position | (300A,012C) | DS | 2C | |

9.1.2.3.9 SOP Common Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------|-------------|----|------|---------|
| SOP Class UID | (0008,0016) | UI | 1 | |
| SOP Instance UID | (0008,0018) | UI | 1 | |

9.1.2.4 RT Beams Treatment Record IOD

| IE | Module | Used |
|------------------|------------------------------------|------|
| Patient | Patient Module | Yes |
| Study | General Study Module | Yes |
| Series | RT Series Module | Yes |
| Equipment | General Equipment Module | No |
| Treatment Record | RT General Treatment Record Module | Yes |
| | RT Patient Setup Module | Yes |
| | RT Treatment Machine Record Module | Yes |
| | RT Beams Session Record Module | Yes |
| | SOP Common Module | Yes |

9.1.2.4.1 Patient Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------------|-------------|----|------|---------|
| Patient's Name | (0010,0010) | PN | 2 | |
| Patient ID | (0010,0020) | LO | 2 | |
| Patient's Birth Date | (0010,0030) | DA | 2 | |
| Patient's Sex | (0010,0040) | CS | 2 | |

9.1.2.4.2 General Study Module

| Attribute name | Tag | Vr | Type | Comment |
|--------------------|-------------|----|------|---------|
| Study Instance UID | (0020,000D) | UI | 1 | |

9.1.2.4.3 RT Series Module

| Attribute name | Tag | Vr | Type | Comment |
|---------------------|-------------|----|------|---------|
| Modality | (0008,0060) | CS | 1 | |
| Series Instance UID | (0020,000E) | UI | 1 | |

9.1.2.4.4 RT General Treatment Record Module

| Attribute name | Tag | Vr | Type | Comment |
|-----------------|-------------|----|------|---------|
| Instance Number | (0020,0013) | IS | 1 | |
| Treatment Date | (3008,0250) | DA | 2 | |
| Treatment Time | (3008,0251) | TM | 2 | |

| | | | | |
|------------------------------|-------------|----|---|--|
| Referenced RT Plan Sequence | {300C,0002} | SQ | 2 | |
| >Referenced SOP Class UID | {0008,1150} | UI | 1 | |
| >Referenced SOP Instance UID | {0008,1155} | UI | 1 | |

9.1.2.4.5 RT Patient Setup Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------------|-------------|----|------|---------|
| Patient Setup Sequence | {300A,0180} | SQ | 1 | |
| >Patient Setup Number | {300A,0182} | IS | 1 | |

9.1.2.4.6 RT Treatment Machine Record Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------------------|-------------|----|------|---------|
| Treatment Machine Sequence | {300A,0206} | SQ | 1 | |

9.1.2.4.7 RT Beams Session Record Module

| Attribute name | Tag | Vr | Type | Comment |
|---|-------------|----|------|------------|
| Primary Dosimeter Unit | {300A,00B3} | CS | 1 | |
| Treatment Session Beam Sequence | {3008,0020} | SQ | 1 | |
| >Referenced Beam Number | {300C,0006} | IS | 3 | |
| >Beam Name | {300A,00C2} | LO | 3 | |
| >Beam Type | {300A,00C4} | CS | 1 | |
| >Radiation Type | {300A,00C6} | CS | 1 | |
| >Beam Limiting Device Leaf Pairs Sequence | {3008,00A0} | SQ | 1 | |
| >Number of Wedges | {300A,00D0} | IS | 1 | |
| >Current Fraction Number | {3008,0022} | IS | 2 | |
| >Treatment Delivery Type | {300A,00CE} | CS | 2 | |
| >Treatment Termination Status | {3008,002A} | CS | 1 | |
| >Specified Primary Meterset | {3008,0032} | DS | 3 | |
| >Specified Secondary Meterset | {3008,0033} | DS | 3 | |
| >Delivered Primary Meterset | {3008,0036} | DS | 3 | |
| >Delivered Secondary Meterset | {3008,0037} | DS | 3 | |
| >Specified Treatment Time | {3008,003A} | DS | 3 | |
| >Delivered Treatment Time | {3008,003B} | DS | 3 | |
| >Number of Control Points | {300A,0110} | IS | 1 | |
| >Control Point Delivery Sequence | {3008,0040} | SQ | 1 | |
| >Treatment Control Point Date | {3008,0024} | DA | 1 | |
| >Treatment Control Point Time | {3008,0025} | TM | 1 | |
| >Delivered Meterset | {3008,0044} | DS | 1 | |
| >Gantry Angle | {300A,011E} | DS | 1C | |
| >Patient Support Angle | {300A,0122} | DS | 1C | |
| >Table Top Pitch Angle | {300A,0140} | FL | 1C | |
| >Table Top Roll Angle | {300A,0144} | FL | 1C | |
| >Table Top Vertical Position | {300A,0128} | DS | 2C | |
| >Table Top Longitudinal Position | {300A,0129} | DS | 2C | |
| >Table Top Lateral Position | {300A,012A} | DS | 2C | |
| >Tomo Registration | {300D,0010} | LO | 3 | TOMO_HA_01 |
| >Tomo Registration Translation | {300D,10B0} | DS | 3 | |

| | | | | |
|------------------------------------|-------------|----|---|----|
| >Tomo Registration Rotation | (300D,10B1) | DS | 3 | |
| >Tomo bug Hi art 3.6.1 | (300D,0010) | LO | 3 | LO |
| >Tomo Registration Translation Bug | (300D,10B0) | DS | 3 | |
| >Tomo Registration Rotation Bug | (300D,10B1) | SH | 3 | |

9.1.2.4.8 SOP Common Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------|-------------|----|------|---------|
| SOP Class UID | (0008,0016) | UI | 1 | |
| SOP Instance UID | (0008,0018) | UI | 1 | |

9.1.2.5 Spatial Registration IOD

| IE | Module | Used |
|----------------------|------------------------------------|------|
| Patient | Patient Module | Yes |
| Study | General Study Module | Yes |
| Series | General Series Module | Yes |
| | Spatial Registration Series Module | No |
| Frame of Reference | Frame of Reference Module | Yes |
| Equipment | General Equipment Module | No |
| Spatial Registration | Spatial Registration Module | Yes |
| | Common Instance Reference Module | No |
| | SOP Common Module | Yes |

9.1.2.5.1 Patient Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------------|-------------|----|------|---------|
| Patient's Name | (0010,0010) | PN | 2 | |
| Patient ID | (0010,0020) | LO | 2 | |
| Patient's Birth Date | (0010,0030) | DA | 2 | |
| Patient's Sex | (0010,0040) | CS | 2 | |

9.1.2.5.2 General Study Module

| Attribute name | Tag | Vr | Type | Comment |
|--------------------|-------------|----|------|---------|
| Study Instance UID | (0020,000D) | UI | 1 | |

9.1.2.5.3 General Series Module

| Attribute name | Tag | Vr | Type | Comment |
|---------------------|-------------|----|------|---------|
| Modality | (0008,0060) | CS | 1 | |
| Series Instance UID | (0020,000E) | UI | 1 | |
| Patient Position | (0018,5100) | CS | 2C | |

9.1.2.5.4 Frame of Reference Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------------|-------------|----|------|---------|
| Frame of Reference UID | (0020,0052) | UI | 1 | |

9.1.2.5.5 Spatial Registration Module

| Attribute name | Tag | Vr | Type | Comment |
|-----------------|-------------|----|------|---------|
| Content Date | (0008,0023) | DA | 1 | |
| Content Time | (0008,0033) | TM | 1 | |
| Instance Number | (0020,0013) | IS | 1 | |

| | | | | |
|--|-------------|----|----|--|
| Content Label | (0070,0080) | CS | 1 | |
| Registration Sequence | (0070,0308) | SQ | 1 | |
| >Frame of Reference UID | (0020,0052) | UI | 1C | |
| >Referenced Image Sequence | (0008,1140) | SQ | 1C | |
| >Referenced SOP Class UID | (0008,1150) | UI | 1 | |
| >Referenced SOP Instance UID | (0008,1155) | UI | 1 | |
| >Matrix Registration Sequence | (0070,0309) | SQ | 1 | |
| >Matrix Sequence | (0070,030A) | SQ | 1 | |
| >Frame of Reference Transformation Matrix | (3006,00C6) | DS | 1 | |
| >Frame of Reference Transformation Matrix Type | (0070,030C) | CS | 1 | |

9.1.2.5.6 SOP Common Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------|-------------|----|------|---------|
| SOP Class UID | (0008,0016) | UI | 1 | |
| SOP Instance UID | (0008,0018) | UI | 1 | |

9.1.2.6 RT Beams Delivery Instruction RETIRED IOD

| IE | Module | Used |
|-----------|--------------------------------------|------|
| Patient | Patient Module | No |
| | Clinical Trial Subject Module | No |
| Study | General Study Module | Yes |
| | Patient Study Module | No |
| | Clinical Trial Study Module | No |
| Series | General Series Module | Yes |
| | Clinical Trial Series Module | No |
| Equipment | General Equipment Module | No |
| Plan | RT Beams Delivery Instruction Module | Yes |
| | Common Instance Reference Module | No |
| | General Reference Module | No |
| | SOP Common Module | Yes |

9.1.2.6.1 General Study Module

| Attribute name | Tag | Vr | Type | Comment |
|--------------------|-------------|----|------|---------|
| Study Instance UID | (0020,000D) | UI | 1 | |

9.1.2.6.2 General Series Module

| Attribute name | Tag | Vr | Type | Comment |
|---------------------|-------------|----|------|---------|
| Modality | (0008,0060) | CS | 1 | |
| Series Instance UID | (0020,000E) | UI | 1 | |

9.1.2.6.3 RT Beams Delivery Instruction Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------------------|-------------|----|------|--|
| Beam Task Sequence | (0074,1020) | SQ | 1 | |
| >Beam Task Type | (0074,1022) | CS | 1 | |
| >Treatment Delivery Type | (300A,00CE) | CS | 1 | Supported values: TREATMENT, CONTINUATION. |
| >Continuation Start Meterset | (0074,0120) | FD | 1C | |

| | | | | |
|----------------------------|-------------|----|----|--|
| >Continuation End Meterset | (0074,0121) | FD | 1C | |
| >Current Fraction Number | (3008,0022) | IS | 1 | |
| >Referenced Beam Number | (300C,0006) | IS | 1 | |
| Omitted Beam Task Sequence | (300C,0111) | SQ | 3 | |
| >Referenced Beam Number | (300C,0006) | IS | 1 | |
| >Reason for Omission | (300C,0112) | CS | 1 | |

9.1.2.6.4 SOP Common Module

| Attribute name | Tag | Vr | Type | Comment |
|------------------|-------------|----|------|---------|
| SOP Class UID | (0008,0016) | UI | 1 | |
| SOP Instance UID | (0008,0018) | UI | 1 | |

9.1.2.7 Unified Procedure Step RETIRED IOD

| IE | Module | Used |
|------------------------|---|------|
| Unified Procedure Step | SOP Common Module | Yes |
| | Unified Procedure Step Relationship Module | Yes |
| | Unified Procedure Step Scheduled Procedure Information Module | Yes |
| | Unified Procedure Step Progress Information Module | Yes |
| | Unified Procedure Step Performed Procedure Information Module | No |
| | Patient Medical Module | No |
| | Transaction Module | Yes |

9.1.2.7.1 SOP Common Module

| Attribute name | Tag | Vr | Type | Comment |
|--------------------------|-------------|----|------|---------|
| SOP Class UID | (0008,0016) | UI | 1 | |
| SOP Instance UID | (0008,0018) | UI | 1 | |
| Timezone Offset From UTC | (0008,0201) | SH | 3 | |

9.1.2.7.2 Unified Procedure Step Relationship Module

| Attribute name | Tag | Vr | Type | Comment |
|----------------|-------------|----|------|---------|
| Patient's Name | (0010,0010) | PN | | |
| Patient ID | (0010,0020) | LO | | |

9.1.2.7.3 Unified Procedure Step Scheduled Procedure Information Module

| Attribute name | Tag | Vr | Type | Comment |
|---|-------------|----|------|---------|
| Scheduled Station Name Code Sequence | (0040,4025) | SQ | | |
| >Code Value | (0008,0100) | SH | 1C | |
| >Coding Scheme Designator | (0008,0102) | SH | 1C | |
| >Code Meaning | (0008,0104) | LO | 1 | |
| Scheduled Procedure Step Start DateTime | (0040,4005) | DT | | |
| Scheduled Workitem Code Sequence | (0040,4018) | SQ | | |
| >Code Value | (0008,0100) | SH | 1C | |
| >Coding Scheme Designator | (0008,0102) | SH | 1C | |
| >Code Meaning | (0008,0104) | LO | 1 | |
| Study Instance UID | (0020,000D) | UI | | |
| Input Availability Flag | (0040,4020) | CS | 1 | |

9.1.2.7.4 Unified Procedure Step Progress Information Module

| Attribute name | Tag | Vr | Type | Comment |
|--|-------------|----|------|---------|
| Procedure Step State | (0074,1000) | CS | | |
| Procedure Step Progress Information Sequence | (0074,1002) | SO | | |
| >Procedure Step Progress | (0074,1004) | DS | | |
| >Procedure Step Progress Description | (0074,1006) | ST | | |

9.1.2.7.5 Transaction Module

| Attribute name | Tag | Vr | Type | Comment |
|-----------------|-------------|----|------|---------|
| Transaction UID | (0008,1195) | UI | 3 | |

9.1.3 Attribute Mapping

Not applicable

9.1.4 Coerced/Modified Fields

Not applicable

9.2 DATA DICTIONARY OF PRIVATE ATTRIBUTES

All used Private Creators are listed in the table below. Usage of Private Attributes are listed in each module specification.

| Attribute name | Tag | VR | VM | Value |
|-----------------------|-------------|----|----|------------|
| Tomo Registration | (300D,0010) | LO | 1 | TOMO_HA_01 |
| Tomo bug Hi art 3.6.1 | (300D,0010) | LO | 1 | LO |

9.3 CODE TERMINOLOGY AND TEMPLATES

Not applicable

9.4 GRAYSCALE IMAGE CONSISTENCY

Not applicable

9.5 STANDARD EXTENDED/SPECIALIZED/PRIVATE SOP CLASSES

9.5.1 Standard extended SOP Class

9.5.2 Specialized SOP Class

Not applicable

9.5.3 Private SOP Class

Not applicable

9.6 PRIVATE TRANSFER SYNTAXES

Not applicable



CONTACT INFORMATION

RaySearch Laboratories AB (publ) - Head office

P.O. Box 3297

SE-103 65 Stockholm, Sweden

Phone: +46 8 510 530 00

Fax: +46 8 510 530 30

Visiting address:

Sveavägen 44

SE-111 34 Stockholm, Sweden

info@raysearchlabs.com

www.raysearchlabs.com

RaySearch Americas

Phone: +1 877 778 3849

RaySearch France

Phone: +33 975 433 632

RaySearch Korea

Phone: +82 10 2230 2046

RaySearch Belgium

Phone: +32 2 213 83 65

RaySearch Germany

Phone: +49 30 89 36 06 90

RaySearch Singapore

Phone: +65 81 28 59 80

RaySearch China

Phone: +86 137 0111 5932

RaySearch Japan

Phone: +81 3 4405 6902

RaySearch UK

Phone: +44 7508 426 563